

MATERIAL SAFETY DATA SHEET

PRODUCT NAME: Prespotter 03

MANUFACTURER:

Acorn Biotechnical Corporation

DATE PREPARED: 03.01.02

P.O. Box 925521
Houston, Texas 77292
(800) 982-1187

HAZARDOUS INGREDIENTS	C.A.S. #	OSHA PEL/TLV	ACGIH TLV	STEL	OTHER LIMITS
Isopropanol	67-63-0	980 mg/m ³	980 mg/m ³	1225 mg/m ³	N/A

HEALTH - 1

FIRE - 1

REACTIVITY - 0

SPECIFIC HAZARD - None

BOILING POINT - > 200°F.

FLASH POINT - N/D

METHOD USED -

SPECIFIC GRAVITY - 1.07

EVAPORATION RATE (WATER = 1) - 1

% VOLATILE - 4

PHYSICAL DESCRIPTION - Slightly turbid mobile liquid with fragrance

FLAMMABLE LIMITS: LEL - N/A

UEL - N/A

EXTINGUISHING MEDIA: Water fog, chemical foam, carbon dioxide

SPECIAL FIRE FIGHTING PROCEDURE: Firefighters should wear a self-contained breathing apparatus as in all chemical fires.

UNUSUAL FIRE AND EXPLOSION HAZARDS: None

UNSTABLE -

STABLE - X

HAZARDOUS DECOMPOSITION PRODUCTS: May produce oxides of carbon and nitrogen.

CONDITION TO AVOID: Temperatures above 110°F.

MATERIALS TO AVOID: Strong acids and bases. Active halogen compounds

HAZARDOUS POLYMERIZATION: Will not occur.

POSSIBLE ROUTES OF ENTRY:

INHALATION - Yes

INGESTION - Yes

SKIN - No

ACUTE HEALTH HAZARD: May produce eye irritation. May produce gastric distress upon ingestion.

CHRONIC HEALTH HAZARD: Repeated exposure may lead to allergic sensitization.

CARCINOGENICITY: NTP - No

IARC MONOGRAPHS - No

OSHA REGULATED - No

ORAL LD/50 - N/D

DERMAL IRRITATION - Slight

EYE IRRITATION - Slight

SYMPTOMS OF OVEREXPOSURE

EYE - Slight irritation

SKIN - Slight irritation will occur if not removed with water or if saturated clothing is not removed and cleaned before wearing.

INHALATION - Irritation of nose and throat. Repeated inhalation of aerosols may produce allergic sensitization.

INGESTION - May produce gastric distress and throat irritation.

CONDITIONS AGGRAVATED BY EXPOSURE - Possible increase in the incidence of allergic asthma

FIRST AID

SKIN: Flush area with a large amount of water.

EYES: Flush eye with water for 15 minutes. If irritation persists, see a doctor.

INHALATION: Avoid producing aerosols. If an aerosol is produced leave area until mist has settled. If chronic irritation is suspected, see a doctor.

INGESTION: Drink large quantity of water or milk. DO NOT INDUCE VOMITING. If distress or irritation persists, see a doctor.

SPILL SAFETY: Absorb material on paper. Clean residue with water.

DISPOSAL INFORMATION: Dispose of paper in appropriate receptacle. Flush diluted material to sewer.

STORAGE AND HANDLING: Store in a cool place.



PERSONAL HYGIENE: Wash all exposed surfaces with soap and water. Do not allow material to contact open wounds. If such contact occurs, wash with soap and water and apply a superficial disinfectant.

PROTECTION - None required TYPE - N/A

VENTILATION:

GENERAL MECHANICAL - X

LOCAL -

SPECIAL -

PROTECTIVE EQUIPMENT REQUIRED:

RUBBER GLOVES -

RUBBER APRON -

RUBBER BOOTS -

CHEMICAL SPLASH GOGGLES -

FACE SHIELD -

SHIPPING INFORMATION:

DOT LABEL - None required

U.N. NUMBER -

SPECIAL INFORMATION:

N/A = NOT APPLICABLE

N/K = NOT KNOWN

N/D = NOT DONE

OTHER COMMENTS

SARA TITLE REPORTING REQUIREMENTS

Section 302,304

Section 311,312

Section 313

CERCLA

Reporting not required

Reporting not required for inventory above TPQ

Reporting may be required for manufacturing use

Same as section 304

The data and recommendations on this document are believed to be accurate. However, no warranty is expressed or implied in regard to these data or the results obtained from the use of the data. In addition Acorn Biotechnical Corporation assumes no responsibility for injury to customers and third parties caused by the material if reasonable safety procedures such as are outlined in this document are not followed. The information in this document is provided by Acorn Biotechnical Corporation as a service to its customers for use in complying with the Occupational Safety and Health Act of 1970. Any other use is prohibited.

CREATED BY SMK

